



**ClariVein.IC**  
Infusion Catheter

**PREPARE. ACCESS. DELIVER.™**

# TOOLS FOR VASCULAR INTERVENTIONS

## PREPARE

Along with standard procedural kits, Merit can create kits to meet your specific clinical requirements. Combining the tools you need into one convenient kit can help increase procedural efficiency and improve inventory management. Our kits combine high quality products for your venous procedures, including safety products, sheaths, and guide wires.




ClariPak™

### KIT CONTENTS:

- Aquasonic® Gel
- Bowls
- Chloraprep™
- Drapes
- Elastic Bandages
- Futura® Safety Scalpel
- Gauze pads
- Gowns
- Hemostat
- Medication Labels
- Needles (18G, 27G)
- Probe Cover Kit
- Skin Marker
- Syringes
- Towels
- Tray

## ACCESS


A small access site minimizes bleeding and trauma to surrounding tissue. Merit micro introducer kits include a coaxial introducer/dilator, a 21-gauge needle, and an 0.018" guide wire. The introducers and dilators provide smooth transitions. A stiffened dilator is also available, to provide additional support during placement in tortuous anatomy. Select kits include a hemostasis valve adapter. All items are available separately.





# DELIVER

ClariVein® IC is an innovative device indicated for the infusion of physician specified agents in the peripheral vasculature. This unique infusion catheter has a 360-degree rotating wire tip at the distal end and does not require thermal energy or tumescence. The device is fully disposable and can be used in an outpatient setting. Using ClariVein IC may offer the patient the benefit of minimal pain and brief recovery time. ^



# TOOLS FOR VASCULAR INTERVENTIONS

## ClariVein®

Catalog Number	Catheter Length
45-018-E2S	45 cm
65-018-E4S	65 cm
85-018-E140S	85 cm

## Procedure Pack

Catalog Number	Description
K22C-00001	ClariPak™ Procedure Pack

Contact your local Merit Sales Representative to design a custom procedural pack for your needs.

## Micro Introducer Kits

Catalog Number	French Size	Length	Guide Wire Diameter/Length	Guide Wire Core/Tip	Needle Gauge/Length	Qty/Box
----------------	-------------	--------	----------------------------	---------------------	---------------------	---------

### STANDARD 7 CM NEEDLE KITS MERIT **MAK**

MAK401	4	10 cm	0.018" (0.46 mm) x 40 cm	SS/SS	21G x 7 cm	10
MAK401NPD	4	10 cm	0.018" (0.46 mm) x 40 cm	NT/PD	21G x 7 cm	10
MAK401N	4	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 7 cm	10
MAK501	5	10 cm	0.018" (0.46 mm) x 40 cm	SS/SS	21G x 7 cm	10
MAK501NPD	5	10 cm	0.018" (0.46 mm) x 40 cm	NT/PD	21G x 7 cm	10
MAK501N	5	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 7 cm	10

### ECHO-ENHANCED 7 CM NEEDLE KITS

MAK401-E	4	10 cm	0.018" (0.46 mm) x 40 cm	SS/SS	21G x 7 cm	10
MAK401NPD-E	4	10 cm	0.018" (0.46 mm) x 40 cm	NT/PD	21G x 7 cm	10
MAK401N-E	4	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 7 cm	10
MAK501-E	5	10 cm	0.018" (0.46 mm) x 40 cm	SS/SS	21G x 7 cm	10
MAK501NPD-E	5	10 cm	0.018" (0.46 mm) x 40 cm	NT/PD	21G x 7 cm	10
MAK501N-E	5	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 7 cm	10

### STIFFENED DILATOR KITS WITH ECHO-ENHANCED 7 CM NEEDLE MERIT **S-MAK**

S-MAK401	4	10 cm	0.018" (0.46 mm) x 40 cm	SS/SS	21G x 7 cm	10
S-MAK401NPD	4	10 cm	0.018" (0.46 mm) x 40 cm	NT/PD	21G x 7 cm	10
S-MAK401N	4	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 7 cm	10
S-MAK501	5	10 cm	0.018" (0.46 mm) x 40 cm	SS/SS	21G x 7 cm	10
S-MAK501NPD	5	10 cm	0.018" (0.46 mm) x 40 cm	NT/PD	21G x 7 cm	10
S-MAK501N	5	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 7 cm	10

### MINI AND STIFFENED ACCESS KITS, INCLUDING HEMOSTASIS VALVE ADAPTER MERIT **PAK** PEDAL ACCESS KIT

PAKM4NPD7	4	10 cm	0.018" (0.46 mm) x 40 cm	NT/PD	21G x 7 cm	10
PAKS4NPD7	4	10 cm	0.018" (0.46 mm) x 40 cm	NT/PD	21G x 7 cm	10
PAKM5NPD7	5	10 cm	0.018" (0.46 mm) x 40 cm	NT/PD	21G x 7 cm	10
PAKS5NPD7	5	10 cm	0.018" (0.46 mm) x 40 cm	NT/PD	21G x 7 cm	10

SS = Stainless Steel NT = Nitinol PD = Palladium PT = Platinum

^ Elias, S., Lam, Y. L., & Wittens, C. H. A. (2013, March). Mechanochemical Ablation (MOCA): Status and Results [greater than two years data]. Phlebology, Vol. 28 (suppl. 1), pp. 10-14. DOI: 10.1177/026835513477787

All trademarks and registered trademarks are the property of their respective owners.

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.253.1600  
1.800.35.MERIT

Merit Medical Europe, Middle  
East, & Africa (EMEA)  
Amerikalaan 42, 6199 AE  
Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733